

# AED/Bleed Kit Inspection Form

Inspected By:	Location of Inspection (Building):
Date of Inspection	AED Serial Number:

Use the following maintenance checklist when you check your unit.

**Check the following Pass or Fail**

<b>Check the following</b>	<b>PASS</b>	<b>FAIL</b>	<b>Comments</b>
Is the green check showing that the unit is ready to use?	<input type="checkbox"/>	<input type="checkbox"/>	
Is the unit clean, undamaged, and free of excessive wear?	<input type="checkbox"/>	<input type="checkbox"/>	
Are there any cracks or loose parts in the housing?	<input type="checkbox"/>	<input type="checkbox"/>	
Verify that electrodes are within their expiration date.	<input type="checkbox"/>	<input type="checkbox"/>	
Verify electrodes are connected to the unit and sealed in their package. Replace if expired.	<input type="checkbox"/>	<input type="checkbox"/>	Exp. Date: